REF FLOW-TIM

PL INSTRUKCJA OBSŁUGI
CS NÁVOD K OBSLUZE
DA BRUGSANVISNING
DE BENUTZERHANDBUCH
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
EN USER MANUAL
ET KASUTUSJUHEND
HR KORSNIČKE UPUTE
HU FELHASZNÁLÓI ÚTMUTATÓ
IT MANUALE D'USO
LT NAUDOJIMO INSTRUKCIJA
LV LIETOTĀJA INSTRUKCIJA
PT MANUAL DO UTILIZADOR
RO MANUAL DE INSTRUCTIUNI
SK NÁVOD NA POUŽITIE
SL NAVODILA ZA UPORABO
SV BRUKSANVISNING

PL INSTRUKCJA OBSŁUGI.....	4
CS NÁVOD K OBSLUZE.....	10
DA BRUGSANVISNING.....	15
DE BENUTZERHANDBUCH.....	20
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	26
EN USER MANUAL	32
ET KASUTUSJUHEND	37
HR KORSNIČKE UPUTE	42
HU FELHASZNÁLÓI ÚTMUTATÓ	47
IT MANUALE D'USO	53
LT NAUDOJIMO INSTRUKCIJA	59
LV LIETOTĀJA INSTRUKCIJA	65
PT MANUAL DO UTILIZADOR.....	70
RO MANUAL DE INSTRUCȚIUNI	76
SK NÁVOD NA POUŽITIE.....	81
SL NAVODILA ZA UPORABO.....	86
SV BRUKSANVISNING.....	91

ELEMENTY WÓZKA

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ELEMENTS OF THE WHEELCHAIR

- | | |
|--|---|
| 1. Handgrip | 9. Solid, rear 22" wheel with handrims* |
| 2. Height adjustment | 10. Footrest lever |
| 3. Backrest folding mechanism | 11. Footrest with calf strap |
| 4. Backrest | 12. Center of gravity adjustment |
| 5. Seat with safety belt | 13. Solid, front 8" wheel |
| 6. Armrest flip-up function | 14. Attendant brake |
| 7. Parking brake (hand brake) | 15. Anti-tip wheel |
| 8. Quick-release function and center of gravity adjustment | |

** available also in version 24" rear wheels (only for seat size 46 cm)*

1. Characteristics

The wheelchair is made of a lightweight aluminum frame. It features flip-up armrests and an optional detachable upholstery (for backrest and seat) made of a breathable mesh material. Thanks to the height adjustment and quick release function on rear wheels and center of gravity adjustment on front and rear wheels, wheelchair can be adjusted to the user's individual needs. It is equipped with foldable backrest mechanism, attendant brakes, safety belt/straps and anti-tip wheels. Wheelchair is available in version on 22" or 24" (only for seat size 46 cm) rear solid wheels and 8" front solid wheels.

The set includes: wheelchair frame with wheels, seat and backrest, safety belt, assembly keys.

Optional accessories: anti-tip wheels, footrests with safety straps, mesh upholstery.

2. Application

The wheelchair is designed for people who have problems with unassisted moving and for patients with dysfunctions and reduced mobility. This self-propelled wheelchair can also be pushed by an accompanying person. Thanks to the ability to adjust the settings to the individual needs of the user, the wheelchair provides better control, greater stability, and a more dynamic, active ride on everyday bases. The wheelchair is suitable for indoor and outdoor use.

3. Assembly

Step 1. Unfolding the wheelchair

Place the wheelchair on the ground. Unfold the wheelchair using the armrests. Then, press the center of the seat so that all the elements fit together.

Step 2. Backrest, seat and height adjustment

Lift the backrest upwards and lock it with backrest locks which are located below the hand grips. Fix the upholstery on the wheelchair. Untight mounting screws under handles, select desired height and tighten them afterwards. Make sure both handles are on the same height level.

Step 3. Footrests

Fix the footrest on the outer side of the frame. Then turn the footrest inwards until it is locked.

Step 4. Armrests

The armrests are tiltable. To raise the armrest, press the button located on the frame (under armrests) of the wheelchair.

Step 5. Height adjustment of footrests

Adjust the footrests using the screw located on the tube. Enough space should be left between the footplates and the ground as footplates set too low may cause difficulties in climbing hills and curbs. Adjust the footplates so that the angle between trunk and thighs is 90° (optimum sitting position for the user).

Step 6. Assembly of safety straps

Attach the safety straps to the footrests to prevent the user's legs from sliding down the footplates.

4. Before use

Step 1. Checking the condition of the wheelchair

Check the condition of the brakes, construction and tires. Check if all the components are assembled and adjusted properly.

Note! Do not use the wheelchair if any defects are noticed.

Step 2. Using parking brakes

Lock the wheels when getting in or out of the wheelchair.

Step 3. Footplates

The footplates should be set in vertical position when getting in or out of the wheelchair. Do not use the footplates as a support while getting in or out of the wheelchair, due to the risk of falling. After moving to the wheelchair, unfold the footplates.

Note! Do not stand on the footplates due to the risk of damage.

Step 4. Getting in the wheelchair

The wheelchair should be located as close to the patient as possible. The wheelchair should touch the bed or other seat being used. If the user cannot change their seat independently, transfer the user to the wheelchair.

Step 5. Position of the user

The user should sit in an upright position and their knees should be bent at the angle of 90° .

Wheelchair is ready to use.

FOLDING

The wheelchair can be re-folded in order to be transported and stored more easily.

- Swing away the footplates and set them in vertical position.
- Detach the upholstery.
- Use the handles to lift the seat upwards until the frame of the wheelchair is folded.

- To fold the backrest, simultaneously press locking elements which are located under the hand grips at the back of the backrest and pull the backrest downwards.

PARKING BRAKE

To lock the wheels, push the lever forward until the wheel is fully blocked (separately for left and right wheel). Please, note that the parking brake only serves to lock the wheels when the wheelchair is at a complete stop.

PARKING BRAKES ADJUSTMENT

The proper functioning of the parking brake depends on its adjustment and position. It should be mounted in adequate distance from the back wheel, so that it is blocked firmly and the lever does not face excessive resistance.

- In case of excessive resistance, move the parking brake towards the front using the locking screws. Set the parking brake in a position that allows it to block the wheels easily and without excessive resistance.
- When the parking brake does not touch or lock the wheel properly, follow the above instruction by moving the handle in the opposite direction.

ATTENDANT BRAKES

The levers of the brakes are located behind the backrest, under both grips.

- In order to stop the wheelchair, squeeze the lever of the brake. To release the brake, let go of the levers. Perform this activity smoothly and with appropriate force.
- The brakes are equipped with parking locks. In order to use them, squeeze it and pull out the locking element.

5. Technical data

Weight 43 (22"): 11,00 kg 46 (22"): 11,47 kg 46 (24"): 11,97 kg 48 (22"): 11,55 kg	Overall width: 61 / 64 / 66 cm
	Seat dimensions: 43 / 46 / 48 x 40 cm
	Overall depth: 106 cm
	Overall height: 100-106 cm
	Seat height from the ground: 50 cm
	Folded width: 106 x 29 x 92 cm
	Maximum capacity: 115 kg

6. Notes

Timago International Group does not bear responsibility of improper usage of the product, failure to comply safety regulations and misuse.

7. Cleaning and maintenance

Carry out the maintenance of the wheelchair components regularly:

- The upholstery and the cushion can be cleaned with damp (not wet) cloth and cold or warm water, without the use of detergents. Difficult stains can be cleaned with a sponge. Do not soak the upholstery for too long. Dry away.
- Plastic elements should be cleaned with damp cloth possibly with the use of appropriate detergents for plastic elements.

8. Storage and transport conditions

Product should not be exposed to direct sunlight, water or moisture.

9. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

10. Medical device incidents

Any serious incident involving a medical device must be reported to the manufacturer and to the competent authority of the European Union member state where the user or patient resides.

11. Labels

	Reference number		Manufacturer
	LOT number		UDI code
	Serial number		Medical device
	Note		Please read the instruction
	The product meets the legal requirements of the European Union for medical devices (MDR 2017/745)		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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